


Microtracer® color card

This color reference
chart has been
provided as a service
to our U.S. customers.

Elanco

Products represented are labeled for use in the United States.




Elanco

Microtracer FS-Red/Natural Yellow (Rumensin®)
Developer: water or 50% alcohol




Phibro

Micotracer RF-Brown-Lake (Terramycin®)
Developer: 7% aqueous sodium carbonate*



Phibro

Microtracer RF-Violet (Red/Blue)-Lake (Stafac®)
Developer: 7% aqueous sodium carbonate*



Zinpro

Microtracer F-Green/Yellow (Zinpro® 120)
Developer: 50% aqueous ethanol



Zinpro

Microtracer F-Red/Yellow-Lake (Zinpro® AvailaMins)
Developer: 7% aqueous sodium carbonate*

* Color spots will not develop using water or water/alcohol mix. Must use 7% aqueous sodium carbonate to develop.

Micro-Tracers Inc. supplied sample results and photos. Colors have been reproduced as accurately as possible.
Some variation may occur.

Elanco™, Rumensin® and the diagonal bar are all trademarks owned or licensed by Eli Lilly and Company, its subsidiaries or affiliates. Microtracer® is a registered trademark of Micro-Tracers Inc. Terramycin® and Stafac® are trademarks of Phibro Animal Health. All other trademarks are the property of their respective owners.

© 2015 Elanco Animal Health. All rights reserved.

USDBURUM00195

Elanco

