

Trigger-Lock Syringe



Attributes

- Allows accurate delivery of vaccines and injectable antibiotics
- Reduces likelihood of unplanned injections
- Trigger-lock mechanism prevents movement of plunger until manually released
- Flexibility in dosing; injection volume can be set in 0.5mL increments from 0.5mL to 5mL
- Constructed of durable plastic



MADE IN AUSTRALIA
N J PHILLIPS Pty Limited ABN 36 000 082 002



DURVET, INC. ANIMAL HEALTH PRODUCTS
P.O. Box 279 · Blue Springs, Missouri 64014 · 800/821-5570 · 816/229-9101
Fax 816/224-3080 · E-mail: info@durvet.com · www.durvet.com